

## HRP-conjugated Goat Anti-Rabbit IgG

Catalog NO.: SA003

For research use only.

Overview

Product name HRP-conjugated Goat Anti-Rabbit IgG

Applications WB,IHC,ELISA

Recommended dilutions WB 1:5000-1:10000

IHC 1:100-1:200

ELISA Use at an assay dependent concentration

NOTE: Optimal dilutions should be determined by the end user.

Size 100µL

**Storage** PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20°C for one year. Avoid repeated freeze-

https://www.elkbiotech.com T:86-27-59760950 E: elkbio@elkbiotech.com

thaw cycles.

Note

This product is for research use only,not for diagnostic or therapeutic use.