



ELK Biotechnology

CY3-conjugated Donkey Anti- Goat IgG

Catalog NO. : SA009

For research use only.

Overview

Product name	CY3-conjugated Donkey Anti- Goat IgG
Recommended dilutions	1:20-1:100 for most applications NOTE: Optimal dilutions should be determined by the end user.
Size	100µL
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Note	This product is for research use only,not for diagnostic or therapeutic use.
Fluorophore	Cy3
Absorption Peak	550nm
Emission Peak	570nm