



ELK Biotechnology

MAKE RESEARCH EASIER

Microsample ELISA kit

2000+ Kit

25 μ L only



Significant Sample Saving

25 μ L of sample required only! (traditional is 100 μ L), reduced by 75% directly. Precious clinical samples, model animal sera, and mouse sera with small blood collection volume will no longer be in short supply, allowing every sample of yours to be fully utilized!

High Sensitivity

pg/mL Level LOD with high-affinity antibodies and an amplified system, enabling clear detection of even trace low-abundance targets. Uncover subtle biological changes and achieve more significant results for high-impact publications!

High Specificity

Cross-reactivity **<5%**, so you can eliminate high background interference.

Outstanding Precision

Intra- and inter-assay **CV < 8.5%**, delivering reliable and consistent results with excellent reproducibility. Trustworthy data provides a solid foundation for your outstanding research outcomes.

ELK Biotechnology CO.,LTD

10410 Corporate Dr, STE12 Sugar Land, TX 77478, USA

www.elkbiotech.com

Research area of ELK Microsample ELISA kits

ELK Biotechnology Microsample ELISA Kits feature an extensive product line, successfully validated for over **800** hot-selling targets across **human, mouse, and rat species**, fully empowering your cross-disciplinary research.



- **Immunology and Inflammation Research**
TNF- α , IL-6, IL-12, CRP, chemokines, IFN- γ , and other inflammatory mediators.
- **Cancer and Tumor**
EGFR, CD19, CEA, PSA, and various growth factors, facilitating research in tumor mechanisms and diagnostics.
- **Metabolic Endocrinology**
GLP1R, PPAR α/δ , LEPR, PTH – enabling in-depth research on diabetes, obesity, and other metabolic diseases.
- **Cardiovascular Research**
Cardiac Troponin (cTnI/T), BNP, PCSK9, NLRP3 – Key tools for precision diagnosis of cardiovascular diseases.
- **Stem Cell and Cell Therapy:** CD34, CD19, GPC3 – delivering critical data for regenerative medicine research.

Cat.NO	Product name	Detection range
ELK1157MS	Mouse IL6(Interleukin 6) Microsample ELISA Kit	7.82-500 pg/mL
ELK1156MS	Human IL6(Interleukin 6) Microsample ELISA Kit	7.82-500 pg/mL
ELK1387MS	Mouse TNFa(Tumor Necrosis Factor Alpha) Microsample ELISA Kit	15.63-1000 pg/mL
ELK1190MS	Human TNFa(Tumor Necrosis Factor Alpha) Microsample ELISA Kit	15.63-1000 pg/mL
ELK10900MS	Mouse TRAF6(TNF Receptor Associated Factor 6) Microsample ELISA Kit	0.32-20 μ g/mL
ELK3775MS	Mouse sIL-6R(SolubleInterleukin 6 Receptor) Microsample ELISA Kit	0.32-20 ng/mL
ELK2507MS	Mouse IL-12p40(Interleukin 12 p40) Microsample ELISA Kit	15.63-1000 pg/mL
ELK9605MS	Mouse CXCL15(Chemokine C-X-C-Motif Ligand 15) Microsample ELISA Kit	78.13-5000 pg/mL
ELK1010MS	Mouse ECF(Eosinophil Chemotactic Factor) Microsample ELISA Kit	7.82-500 pg/mL
ELK10778MS	Mouse LECT2(Leukocyte Cell Derived Chemotaxin 2) Microsample ELISA Kit	1.57-100 ng/mL
ELK10579MS	Mouse SLC(Secondary Lymphoid Tissue Chemokine) Microsample ELISA Kit	31.25-2000 pg/mL
ELK2237MS	Mouse CD19(Cluster Of Differentiation 19) Microsample ELISA Kit	0.16-10 ng/mL
ELK2421MS	Mouse EGFR(Epidermal Growth Factor Receptor) Microsample ELISA Kit	93.75-6000 pg/mL
ELK9614MS	Mouse PSA(ProstateSpecific Antigen) Microsample ELISA Kit	0.16-10 ng/mL
ELK6119MS	Mouse GLP1R(Glucagon Like Peptide 1 Receptor) Microsample ELISA Kit	0.16-10 ng/mL
ELK2009MS	Mouse PPARa(Peroxisome Proliferator Activated Receptor Alpha) Microsample ELISA Kit	0.32-20 ng/mL
ELK9212MS	Rat hs-cTnI(High Sensitivity Cardiac Troponin I) Microsample ELISA Kit	15.63-1000 pg/mL
ELK6207MS	Mouse cTnT/TNNT2(Troponin T Type 2, Cardiac) Microsample ELISA Kit	15.63-1000 pg/mL
ELK2136MS	Mouse CTNNb1(Catenin Beta 1) Microsample ELISA Kit	0.16-10 ng/mL
ELK2237MS	Mouse CD19(Cluster Of Differentiation 19) Microsample ELISA Kit	0.16-10 ng/mL
ELK0918MS	Mouse GPC3(Glypican 3) Microsample ELISA Kit	0.16-10 ng/mL
ELK10393MS	Mouse LEPR(Leptin Receptor) Microsample ELISA Kit	3.13-200 ng/mL