


ELK Biotechnology

Cytokeratin 6 Mouse mAb Catalog NO.: EM1080 For research use only.

Overview

Product name	Cytokeratin 6 Mouse Monoclonal antibody
Source	Mouse
Applications	WB IHC
Species reactivity	Human
Recommended dilutions	WesternBlot:1/1000 Immunohistochemistry:1/200 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1 mg/ml
Observed band	56-60kDa
GeneID (Human)	286887
Human Swiss-Prot No.	P02538
Cellular localization	Intermediate filament
Alternative Names	N/A
Background	Cytokeratins are proteins of keratin-containing intermediate filaments found in the intracytoplasmic cytoskeleton of epithelial tissue.



Western blot analysis of) Hela 2) ${\rm HepG2}$ with Cytokeratin Mouse mAb diluted at:2000.