


ELK Biotechnology P53(Mono Methyl Lys370) Mouse mAb Catalog NO.: EM1198 For research use only.

Overview

Product name	P53(Mono Methyl Lys370) Mouse Monoclonal antibody
Source	Mouse
Applications	IHC
Species reactivity	Human Rat Mouse
Recommended dilutions	Immunohistochemistry:1/500-1000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
lsotype	lgG1
Clonality	Monoclonal
Concentration	1 mg/ml
Concentration Observed band	1 mg/ml 53kDa
	•
Observed band	53kDa
Observed band GenelD (Human)	53kDa 7157
Observed band GenelD (Human) Human Swiss-Prot No.	53kDa 7157 P04637



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using P53(Mono Methyl Lys370) (EM1198) Mouse mAb diluted at:500.