




# ELK Biotechnology

RFP-Tag Mouse pAb  
Catalog NO.: EA010  
For research use only.

## Overview

---

Product name	RFP-Tag Mouse polyclonal antibody
Source	Mouse
Applications	WB, IP, IF
Species reactivity	N/A
Recommended dilutions	<b>WesternBlot:1/2000-5000</b> <b>Immunohistochemistry:1/200</b> <b>Immunoprecipitation:1/100</b> <b>NOTE: Optimal dilutions should be determined by the end user.</b>
Immunogen	Recombinant Protein
Species	N/A
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	N/A
GeneID (Human)	N/A
Human Swiss-Prot No.	N/A
Cellular localization	N/A
Alternative Names	DsRed, Red fluorescent protein, Red fluorescent protein drFP583
Background	The red fluorescent protein cloned from Discosoma coral (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein(GFP) from the Aequorea jellyfish.



Western blot analysis of RFP Recombinant protein with EA010 diluted at  
1、 1:5,000  
2、 1:10,000