

BFP-Tag Mouse pAb Catalog NO.: EA032 For research use only.

Overview

Product name BFP-Tag Mouse polyclonal antibody

Source Mouse

Applications WB

Species reactivity N/A

Recommended dilutions WesternBlot:1/2000-5000

Immunohistochemistry:1/200 Immunofluorescence:1/200 Immunoprecipitation:1/200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Recombinant Protein

Species N/A

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

Isotype IgG

Clonality Polyclonal

Concentration 1 mg/ml

Observed band N/A

GenelD (Human) N/A

Human Swiss-Prot No. N/A

Cellular localization N/A

Alternative Names BFP, EBFP, Blue Fluorescent Protein, Enhanced Blue Fluorescent Protein

Background Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue

fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.



Western blot analysis of EBFP recombinant protein with EA032 diluted at 1:5,000.