



ELK Biotechnology

mStrawberry-Tag Mouse pAb

Catalog NO.: EA033

For research use only.

Overview

Product name	mStrawberry-Tag Mouse polyclonal antibody
Source	Mouse
Applications	WB
Species reactivity	N/A
Recommended dilutions	WesternBlot:1/5000 Immunohistochemistry:1/200 Immunofluorescence:1/200 Immunoprecipitation:1/200 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	N/A
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	N/A
GeneID (Human)	N/A
Human Swiss-Prot No.	N/A
Cellular localization	N/A
Alternative Names	mStrawberry
Background	mStrawberry is a bright, monomeric red fluorescent protein which was developed in Dr. Roger Tsien's lab by directed mutagenesis of mRFP.



Western blot analysis of mStawberry recombinant protein with EA033 diluted at 1:5,000.