



ELK Biotechnology

GAL4 Activation Domain Rabbit pAb

Catalog NO.: EA201

For research use only.

Overview

Product name	GAL4 Activation Domain Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	N/A
Recommended dilutions	WesternBlot:1/10000-50000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	N/A
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	N/A
GeneID (Human)	N/A
Human Swiss-Prot No.	N/A
Cellular localization	N/A
Alternative Names	N/A
Background	GAL4 is a transcriptional activator which controls the expression of genes involved in the catabolism of galactose.



Western blot analysis of Recombinant AD protein with GAL4 Activation
Domian Rabbit pAb EA201 diluted at 1:20,000.