

# STRN3 rabbit pAb

Cat No.:ES10750

For research use only

## Overview

Product Name	STRN3 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein AA range: 146-196
Specificity	STRN3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Striatin-3 (Cell cycle autoantigen SG2NA) (S/G2 antigen)
Gene Name	STRN3 GS2NA SG2NA
Cellular localization	Cytoplasm . Membrane ; Peripheral membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	87kD
Human Gene ID	29966
Human Swiss-Prot Number	Q13033
Alternative Names	
Background	alternative products:Additional isoforms seem to exist,caution:Was originally (PubMed:7864889) thought to be nuclear.,function:Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.,similarity:Belongs to the WD repeat striatin family.,similarity:Contains 6 WD repeats.,subunit:Interacts with protein





**ELK Biotechnology**

phosphatase 2A (PP2A) .,



+86-27-59760950

[ELKbio@ELKbiotech.com](mailto:ELKbio@ELKbiotech.com)

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C