



PI16 rabbit pAb

Cat No.:ES10979

For research use only

Overview

Product Name	PI16 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PI16 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Peptidase inhibitor 16 (PI-16) (Cysteine-rich secretory protein 9) (CRISP-9) (PSP94-binding protein)
Gene Name	PI16 CRISP9 PSPBP PSEC0164 UNQ289/PRO328
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	50kD
Human Gene ID	221476
Human Swiss-Prot Number	Q6UXB8
Alternative Names	
Background	function:Putative serine protease inhibitor.,miscellaneous:May serve as a marker following prostatectomy for prostate cancer.,PTM:N-glycosylated.,similarity:Belongs to the CRISP family.,subunit:Interacts with PSP94/MSMB.,tissue specificity:Expressed in prostate, testis, ovary and intestine. Concentrates in





ELK Biotechnology

prostate cancer patient's sera.,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C