



Thymidine kinase Polyclona Antibody

Cat No.:ES12721

For research use only

Overview

Product Name	Thymidine kinase Polyclona Antibody
Host species	Rabbit
Applications	ELISA
Species Cross-Reactivity	Human virus
Recommended dilutions	ELISA(peptide) 1:5000-20000
Immunogen	Synthesized peptide derived from human Thymidine kinase
Specificity	This antibody detects endogenous levels of human Thymidine kinase
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Thymidine kinase
Gene Name	TK
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q9QNF7
Alternative Names	
Background	

