

Raptor (phospho-Ser792) rabbit pAb

Cat No.:ES13461

For research use only

Overview

Product Name	Raptor (phospho-Ser792) rabbit pAb
Host species	Rabbit
Applications	WB;ELISA;IHC
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Immunogen	Synthesized phosho peptide around human Raptor
	(Ser792)
Specificity	This antibody detects endogenous levels of
	Human Mouse Rat Raptor (phospho-Ser792)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Raptor (Ser792)
Gene Name	RPTOR KIAA1303 RAPTOR
Cellular localization	Cytoplasm. Lysosome. Cytoplasmic granule .
	Targeting to lysosomes depends on amino acid
	availability. In arsenite-stressed cells, accumulates in
	stress granules when associated with SPAG5 and
	association with lysosomes is drastically decreased.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	146kD
Human Gene ID	57521
Human Swiss-Prot Number	Q8N122
Alternative Names	Regulatory-associated protein of mTOR (Raptor)
	(p150 target of rapamycin (TOR)-scaffold protein)
Background	This gene encodes a component of a signaling
	pathway that regulates cell growth in response to
	nutrient and insulin levels. The encoded protein
	forms a stoichiometric complex with the mTOR



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



kinase, and also associates with eukaryotic initiation factor 4E-binding protein-1 and ribosomal protein S6 kinase. The protein positively regulates the downstream effector ribosomal protein S6 kinase, and negatively regulates the mTOR kinase. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009],

Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C