

## CN045 rabbit pAb

Cat No.: ES17359

For research use only

## Overview

Product Name CN045 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human CN045 AA

range: 124-174

**Specificity** This antibody detects endogenous levels of CN045

at Human

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name CN045 Gene Name C14orf45

**Cellular localization** Cytoplasm, cytoskeleton, cilium basal body.

Localizes to a polar structure adjacent to the basal

body...

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 80127 Human Swiss-Prot Number Q8ND07


**Alternative Names** 

+86-27-59760950

Background







Western blot analysis of lysates from PC-12 cells, primary antibody was diluted at 1:1000, 4° over night

