

# CJ128 rabbit pAb

Cat No.:ES17440

For research use only

## Overview

Product Name	CJ128 rabbit pAb
Host species	Rabbit
Applications	IHC;IF
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	IHC-p 1: 50-200
Immunogen	Synthesized peptide derived from human CJ128 AA range: 35-85
Specificity	This antibody detects endogenous levels of CJ128 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	CJ128
Gene Name	C10orf128
Cellular localization	Membrane; Single-pass type I membrane protein.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	170371
Human Swiss-Prot Number	Q5T292
Alternative Names	
Background	

