

CHADL rabbit pAb

Cat No.:ES17488

For research use only

Overview

Product Name	CHADL rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000
Immunogen	Synthesized peptide derived from human CHADL AA
	range: 493-543
Specificity	This antibody detects endogenous levels of CHADL at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
l'officiation	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	CHADL
Gene Name	CHADL SLRR4B
Cellular localization	Secreted . Secreted, extracellular space, extracellular
	matrix .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	85kD
Human Gene ID	150356
Human Swiss-Prot Number	Q6NUI6
Alternative Names	Chondroadherin-like protein
Background	




ELKbio@ELKbiotech.com

www.elkbiotech.com

www.eikbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



15 --10 -- Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C