

## CDR2 rabbit pAb

Cat No.: ES17549

For research use only

## Overview

Product Name CDR2 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human CDR2 AA

range: 278-328

**Specificity** This antibody detects endogenous levels of CDR2 at

Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name CDR2

Gene Name CDR2 PCD17
Cellular localization cytoplasm,

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 1039 Human Swiss-Prot Number Q01850


**Alternative Names** 

+86-27-59760950

**Background** 







Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

