



CC037 rabbit pAb

Cat No.:ES17811

For research use only


Overview

Product Name	CC037 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1:500-2000
Immunogen	Synthesized peptide derived from human CC037 AA range: 173-223
Specificity	This antibody detects endogenous levels of CC037 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	CC037
Gene Name	C3orf37 DC12
Cellular localization	Chromosome . Recruited to chromatin following DNA damage (PubMed:30554877). Localizes to replication forks (PubMed:30554877). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	38kD
Human Gene ID	56941
Human Swiss-Prot Number	Q96FZ2
Alternative Names	UPF0361 protein C3orf37
Background	





ELK Biotechnology



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.