

CA167 rabbit pAb


Cat No.:ES17908

For research use only

Overview

Product Name	CA167 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human CA167 AA range: 1015-1065
Specificity	This antibody detects endogenous levels of CA167 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	CA167
Gene Name	C1orf167
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	100506310
Human Swiss-Prot Number	Q5SNV9
Alternative Names	
Background	





Western blot analysis of lysates from VEC cells, primary antibody was diluted at 1:1000, 4° over night

