



BACH1/BRIP1 (phospho-Thr1133) rabbit pAb

Cat No.:ES18105

For research use only

Overview

Product Name	BACH1/BRIP1 (phospho-Thr1133) rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized phospho peptide around human BACH1 and BRIP1 (Thr1133)
Specificity	This antibody detects endogenous levels of Human BACH1 or BRIP1 (phospho-Thr1133)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	BACH1/BRIP1 (Thr1133)
Gene Name	BRIP1 BACH1 FANCI
Cellular localization	Nucleus . Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	135kD
Human Gene ID	83990
Human Swiss-Prot Number	Q9BX63
Alternative Names	Fanconi anemia group J protein (Protein FANCI) (EC 3.6.4.13) (ATP-dependent RNA helicase BRIP1) (BRCA1-associated C-terminal helicase 1) (BRCA1-interacting protein C-terminal helicase 1) (BRCA1-interacting protein 1)

Background

