

ARHGAP42 (phospho-Tyr376) rabbit pAb

Cat No.:ES18261

For research use only

Overview

Product Name	ARHGAP42 (phospho-Tyr376) rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized phospho peptide around human ARHGAP42 (Tyr376)
Specificity	This antibody detects endogenous levels of Human Mouse ARHGAP42 (phospho-Tyr376)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ARHGAP42 (Tyr376)
Gene Name	ARHGAP42
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	95kD
Human Gene ID	143872
Human Swiss-Prot Number	A6NI28
Alternative Names	Rho GTPase-activating protein 42 (Rho GTPase-activating protein 10-like) (Rho-type GTPase-activating protein 42)

Background

