



ALKBH5 rabbit pAb


Cat No.:ES18390

For research use only


Overview

Product Name	ALKBH5 rabbit pAb
Host species	Rabbit
Applications	IHC;IF;WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	IHC-p 1:50-200, WB 1:500-2000
Immunogen	Synthesized peptide derived from human ALKBH5
Specificity	This antibody detects endogenous levels of human ALKBH5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ALKBH5
Gene Name	ALKBH5 ABH5 OFOXD1
Cellular localization	Nucleus speckle .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	51kD
Human Gene ID	54890
Human Swiss-Prot Number	Q6P6C2
Alternative Names	Probable alpha-ketoglutarate-dependent dioxygenase ABH5 (EC 1.14.11.-;Alkylated DNA repair protein alkB homolog 5)
Background	similarity:Belongs to the alkB family.,





Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

