



AKND1 rabbit pAb

Cat No.:ES18422

For research use only


Overview

Product Name	AKND1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human AKND1 AA range: 373-423
Specificity	This antibody detects endogenous levels of AKND1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	AKND1
Gene Name	AKNAD1 C1orf62
Cellular localization	cytoplasm,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	254268
Human Swiss-Prot Number	Q5T1N1
Alternative Names	
Background	This gene encodes a protein which contains a domain found in an AT-hook-containing transcription factor. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2012],





ELK Biotechnology



Western blot analysis of lysates from HuvEc cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.