

AF1L1 rabbit pAb

Cat No.:ES18445

For research use only

Overview


Product Name	AF1L1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1:500-2000
Immunogen	Synthesized peptide derived from human AF1L1 AA range: 180-230
Specificity	This antibody detects endogenous levels of AF1L1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	AF1L1
Gene Name	AFAP1L1
Cellular localization	Cytoplasm . Cell projection, podosome . Cell projection, invadopodium . Cytoplasm, cytoskeleton, stress fiber .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	85kD
Human Gene ID	134265
Human Swiss-Prot Number	Q8TED9
Alternative Names	Actin filament-associated protein 1-like 1 (AFAP1-like protein 1)

Background





ELK Biotechnology



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.