

## ACTY rabbit pAb

## Cat No.:ES18465

For research use only

## Overview

Product Name	ACTY rabbit pAb	
Host species	Rabbit	
Applications	WB	
Species Cross-Reactivity	Human; Mouse	
Recommended dilutions	WB 1: 500-2000	
Immunogen	Synthesized peptide derived from human ACTY AA range: 29-79	
Specificity	This antibody detects endogenous levels of ACTY at Human/Mouse	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	ACTY	
Gene Name	ACTR1B CTRN2	
Cellular localization	Cytoplasm, cytoskeleton. Cytoplasm, cytoskeleton,	
	microtubule organizing center, centrosome.	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band		
Human Gene ID	10120	
Human Swiss-Prot Number	P42025	
Alternative Names		
Background	This gene encodes a 42.3 kD subunit of dynactin, a macromolecular complex consisting of 10 subunits ranging in size from 22 to 150 kD. Dynactin binds to both microtubules and cytoplasmic dynein and is involved in a diverse array of cellular functions, including ER-to-Golgi transport, the centripetal movement of lysosomes and endosomes, spindle	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



positioning, and axonogenesis. This subunit, like ACTR1A, is an actin-related protein. These two proteins, which are of equal length and share 90% amino acid identity, are present in a constant ratio of approximately 1:15 in the dynactin complex. [provided by RefSeq, Aug 2008],

Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4° over night





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C