



ACOT6 rabbit pAb

Cat No.:ES18493

For research use only


Overview

Product Name	ACOT6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human ACOT6 AA range: 116-166
Specificity	This antibody detects endogenous levels of ACOT6 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ACOT6
Gene Name	ACOT6 C14orf42
Cellular localization	[Isoform 1]: Peroxisome . Localization to the peroxisome is uncertain since the potential C-terminal peroxisome targeting signal found in the mouse ortholog is not perfectly conserved. . ; [Isoform 2]: Cytoplasm . Recombinant N-terminally GFP-tagged protei
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	641372
Human Swiss-Prot Number	Q3I5F7
Alternative Names	
Background	





ELK Biotechnology



Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.