

## ABT1 rabbit pAb

## Cat No.:ES18507

For research use only

## **Overview**

ABT1 rabbit pAb
Rabbit
WB
Human; Mouse;Rat
WB 1: 500-2000
Synthesized peptide derived from human ABT1 AA range: 140-190
This antibody detects endogenous levels of ABT1 at Human/Mouse/Rat
Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Store at -20°C. Avoid repeated freeze-thaw cycles.
ABT1
ABT1
Nucleus . Nucleus, nucleolus .
The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Polyclonal
1 mg/ml
29777
Q9ULW3
Basal transcription of genes by RNA polymerase II requires the interaction of TATA-binding protein (TBP) with the core region of class II promoters. Studies in mouse suggest that the protein encoded by this gene likely activates basal transcription from class II promoters by interaction with TBP and the class II promoter DNA. [provided by RefSeq, Jul 2008],



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C.





Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C