

## ABCB9 rabbit pAb

## Cat No.:ES18513

For research use only

## Overview

Product Name	ABCB9 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1:500-2000
Immunogen	Synthesized peptide derived from human ABCB9 AA range: 224-274
Specificity	This antibody detects endogenous levels of ABCB9 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ABCB9
Gene Name	ABCB9 KIAA1520
Cellular localization	Lysosome membrane ; Multi-pass membrane protein . May be located in membrane rafts. Takes an intracellular route from the endoplasmic reticulum (ER), via Golgi and early endosomes to late endosomal and lysosomal compartments (PubMed:30877195)
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	85kD
Human Gene ID	23457
Human Swiss-Prot Number	Q9NP78
Alternative Names	ATP-binding cassette sub-family B member 9 (ATP-binding cassette transporter 9) (ABC transporter 9 protein) (hABCB9) (TAP-like protein) (TAPL)
Background	The membrane-associated protein encoded by this
a	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance as well as antigen presentation. This family member functions in the translocation of peptides from the cytosol into the lysosomal lumen. Alternative splicing of this gene results in distinct isoforms which are likely to have different substrate specificities. [provided by RefSeq, Jul 2011],

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C