



53BP1 (phospho-Ser1618) rabbit pAb

Cat No.:ES18531

For research use only

Overview

Product Name	53BP1 (phospho-Ser1618) rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized phospho peptide around human 53BP1 (Ser1618)
Specificity	This antibody detects endogenous levels of Human Rat 53BP1 (phospho-Ser1618)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	53BP1 (Ser1618)
Gene Name	TP53BP1
Cellular localization	Nucleus . Chromosome . Chromosome, centromere, kinetochore . Localizes to the nucleus in absence of DNA damage (PubMed:28241136). Following DNA damage, recruited to sites of DNA damage, such as double strand breaks (DSBs): recognizes and binds histone H2A monoubiquitinated at 'Lys-15' (H2AK15Ub) and histone H4 dimethylated at 'Lys-20' (H4K20me2), two histone marks that are present at DSBs sites (PubMed:23333306, PubMed:23760478, PubMed:24703952, PubMed:28241136, PubMed:17190600). Associated with kinetochores during mitosis (By similarity). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	213kD
Human Gene ID	7158





ELK Biotechnology

Human Swiss-Prot Number Q12888

Alternative Names Tumor suppressor p53-binding protein 1 (53BP1)
(p53-binding protein 1) (p53BP1)

Background



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C