



DJC24 rabbit pAb

Cat No.:ES18866

For research use only

Overview

Product Name	DJC24 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DJC24 AA range: 39-89
Specificity	This antibody detects endogenous levels of human DJC24
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	DnaJ homolog subfamily C member 24 (CSL-type zinc finger-containing protein 3) (DPH4 homolog)
Gene Name	DNAJC24 DPH4 ZCSL3
Cellular localization	Cytoplasm, cytoskeleton .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	120526
Human Swiss-Prot Number	Q6P3W2
Alternative Names	
Background	

