

NUD15 rabbit pAb

Cat No.:ES18880

For research use only

Overview

Product Name	NUD15 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NUD15 AA range: 62-112
Specificity	This antibody detects endogenous levels of human NUD15
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Probable 8-oxo-dGTP diphosphatase NUDT15 (8-oxo-dGTPase NUDT15) (EC 3.6.1.55) (7,8-dihydro-8-oxoguanine-triphosphatase NUDT15) (MutT homolog 2) (MTH2) (Nucleoside diphosphate-linked moiety X motif 15)
Gene Name	NUDT15 MTH2
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55270
Human Swiss-Prot Number	Q9NV35
Alternative Names	
Background	

