



TRUB1 rabbit pAb

Cat No.:ES18898

For research use only

Overview

| | |
|---------------------------------|---|
| Product Name | TRUB1 rabbit pAb |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human;Mouse;Rat |
| Recommended dilutions | WB 1:1000-2000 |
| Immunogen | Synthesized peptide derived from human TRUB1 AA range: 116-166 |
| Specificity | This antibody detects endogenous levels of human TRUB1 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C . Avoid repeated freeze-thaw cycles. |
| Protein Name | Probable tRNA pseudouridine synthase 1 (EC 5.4.99.-) |
| Gene Name | TRUB1 PUS4 |
| Cellular localization | Nucleus . Cytoplasm, cytosol . Catalyzes pseudouridylation of mRNAs in the nucleus. . |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | 142940 |
| Human Swiss-Prot Number | Q8WWH5 |
| Alternative Names | |
| Background | |

