

ZBED4 rabbit pAb

Cat No.:ES18902

For research use only

Overview

Product Name	ZBED4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZBED4 AA range: 10-60
Specificity	This antibody detects endogenous levels of human ZBED4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger BED domain-containing protein 4
Gene Name	ZBED4 KIAA0637
Cellular localization	Nucleus . Cytoplasm . Photoreceptor inner segment .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9889
Human Swiss-Prot Number	O75132
Alternative Names	
Background	

