

ZDHC4 rabbit pAb

Cat No.:ES18908

For research use only

Overview

Product Name	ZDHC4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZDHC4 AA range: 109-159
Specificity	This antibody detects endogenous levels of human ZDHC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Probable palmitoyltransferase ZDHHC4 (EC 2.3.1.-) (Zinc finger DHHC domain-containing protein 4) (DHHC-4) (Zinc finger protein 374)
Gene Name	ZDHHC4 ZNF374 DC1 UNQ5787/PRO19576
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Golgi apparatus membrane ; Multi-pass membrane protein . Cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55146
Human Swiss-Prot Number	Q9NPG8
Alternative Names	
Background	

