

# CP059 rabbit pAb

Cat No.:ES18969

For research use only

## Overview

|                          |   |
|--------------------------|---|
| Product Name             | CP059 rabbit pAb  |
| Host species             | Rabbit  |
| Applications             | WB  |
| Species Cross-Reactivity | Human;Mouse   |
| Recommended dilutions    | WB 1:1000-2000  |
| Immunogen                | Synthesized peptide derived from human CP059 AA range: 259-309  |
| Specificity              | This antibody detects endogenous levels of human CP059  |
| Formulation              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| Storage                  | Store at -20°C . Avoid repeated freeze-thaw cycles.   |
| Protein Name             | Uncharacterized protein C16orf59  |
| Gene Name                | C16orf59  |
| Cellular localization    | Cell projection, cilium . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole .             |
| Purification             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality                | Polyclonal  |
| Concentration            | 1 mg/ml   |
| Observed band            |   |
| Human Gene ID            | 80178   |
| Human Swiss-Prot Number  | Q7L2K0  |
| Alternative Names        |   |
| Background               |   |

