



F110C rabbit pAb

Cat No.:ES19032

For research use only

Overview

Product Name	F110C rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human F110C AA range: 109-159
Specificity	This antibody detects endogenous levels of human F110C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein FAM110C
Gene Name	FAM110C
Cellular localization	Cytoplasm, cytoskeleton . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle pole . Nucleus . Colocalizes with microtubules during interphase (PubMed:17499476). Detected at the mitotic spindle poles (PubMed:17499476). Colocalizes with AKT1 at the cell cortex (PubMed:19698782). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	642273
Human Swiss-Prot Number	Q1W6H9
Alternative Names	
Background	

