

# PRAM rabbit pAb

Cat No.:ES19148

For research use only

## Overview

Product Name	PRAM rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PRAM AA range: 310-360
Specificity	This antibody detects endogenous levels of human PRAM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	PML-RARA-regulated adapter molecule 1 (PRAM) (PRAM-1)
Gene Name	PRAM1
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84106
Human Swiss-Prot Number	Q96QH2
Alternative Names	
Background	

