

TPIS rabbit pAb

Cat No.:ES19153

For research use only

Overview

Product Name	TPIS rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TPIS AA range: 76-126
Specificity	This antibody detects endogenous levels of human TPIS
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Triosephosphate isomerase (TIM) (EC 5.3.1.1) (Triose-phosphate isomerase)
Gene Name	TPI1 TPI
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7167
Human Swiss-Prot Number	P60174
Alternative Names	
Background	

