



ELP6 rabbit pAb

Cat No.:ES19377

For research use only

Overview

Product Name	ELP6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ELP6 AA range: 22-72
Specificity	This antibody detects endogenous levels of human ELP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Elongator complex protein 6 (Angiotonin-transactivated protein 1) (Protein TMEM103)
Gene Name	ELP6 ATP1 C3orf75 TMEM103
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	54859
Human Swiss-Prot Number	Q0PNE2
Alternative Names	
Background	

