

SCAM4 rabbit pAb

Cat No.:ES19541

For research use only

Overview

Product Name	SCAM4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SCAM4 AA range: 4-54
Specificity	This antibody detects endogenous levels of human SCAM4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Secretory carrier-associated membrane protein 4 (Secretory carrier membrane protein 4)
Gene Name	SCAMP4
Cellular localization	Membrane; Multi-pass membrane protein.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	113178
Human Swiss-Prot Number	Q969E2
Alternative Names	
Background	

