



S41A3 rabbit pAb

Cat No.:ES19782

For research use only

Overview

Product Name	S41A3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human S41A3 AA range: 9-59
Specificity	This antibody detects endogenous levels of human S41A3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Solute carrier family 41 member 3
Gene Name	SLC41A3
Cellular localization	Cell membrane; Multi-pass membrane protein.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	54946
Human Swiss-Prot Number	Q96GZ6
Alternative Names	
Background	

