



CTGE6 rabbit pAb

Cat No.:ES19785

For research use only

Overview

Product Name	CTGE6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CTGE6 AA range: 169-219
Specificity	This antibody detects endogenous levels of human CTGE6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Putative protein cTAGE-6 (cTAGE family member 6)
Gene Name	CTAGE6P CTAGE6
Cellular localization	Membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	340307
Human Swiss-Prot Number	Q86UF2
Alternative Names	
Background	

