



SRCA rabbit pAb

Cat No.:ES19791

For research use only

Overview

Product Name	SRCA rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SRCA AA range: 544-594
Specificity	This antibody detects endogenous levels of human SRCA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Sarcalumenin
Gene Name	SRL
Cellular localization	Sarcoplasmic reticulum lumen . Sarcoplasmic reticulum membrane ; Peripheral membrane protein . May associate with the sarcoplasmic reticulum membrane, via a calcium-dependent mechanism. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	6345
Human Swiss-Prot Number	Q86TD4
Alternative Names	
Background	

