

SERF1 rabbit pAb

Cat No.:ES19800

For research use only

Overview

Product Name	SERF1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SERF1 AA range: 10-60
Specificity	This antibody detects endogenous levels of human SERF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Small EDRK-rich factor 1 (Protein 4F5) (h4F5) (SMA modifier 1)
Gene Name	SERF1A FAM2A SERF1 SMAM1; SERF1B FAM2B SERF1 SMAM1
Cellular localization	Cytoplasm, cytosol . Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	728492
Human Swiss-Prot Number	O75920
Alternative Names	
Background	

