



PDE6C rabbit pAb

Cat No.:ES19801

For research use only

Overview

Product Name	PDE6C rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PDE6C AA range: 474-524
Specificity	This antibody detects endogenous levels of human PDE6C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Cone cGMP-specific 3',5'-cyclic phosphodiesterase subunit alpha' (EC 3.1.4.35) (cGMP phosphodiesterase 6C)
Gene Name	PDE6C PDEA2
Cellular localization	Cell membrane ; Lipid-anchor ; Cytoplasmic side .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	5146
Human Swiss-Prot Number	P51160
Alternative Names	
Background	

