



RIR2 rabbit pAb

Cat No.:ES19802

For research use only

Overview

Product Name	RIR2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RIR2 AA range: 260-310
Specificity	This antibody detects endogenous levels of human RIR2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Ribonucleoside-diphosphate reductase subunit M2 (EC 1.17.4.1) (Ribonucleotide reductase small chain) (Ribonucleotide reductase small subunit)
Gene Name	RRM2 RR2
Cellular localization	Cytoplasm . Nucleus . Localized to the cytoplasm in S phase cells. May localize to the nucleus in G2 phase cells. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	6241
Human Swiss-Prot Number	P31350
Alternative Names	
Background	

