



NARG2 rabbit pAb

Cat No.:ES19804

For research use only

Overview

Product Name	NARG2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NARG2 AA range: 474-524
Specificity	This antibody detects endogenous levels of human NARG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	NMDA receptor-regulated protein 2
Gene Name	NARG2 BRCC1 UNQ3101/PRO10100
Cellular localization	Nucleus . Colocalizes with COIL in subnuclear Cajal and histone locus bodies. Translocates in the LEC complex to Cajal and histone locus bodies at snRNA genes in a ICE1-dependent manner. Associates to transcriptionally active chromatin at snRNA genes.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79664
Human Swiss-Prot Number	Q659A1
Alternative Names	
Background	

