



MSTP9 rabbit pAb

Cat No.:ES19806

For research use only

Overview

Product Name	MSTP9 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MSTP9 AA range: 154-204
Specificity	This antibody detects endogenous levels of human MSTP9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Putative macrophage-stimulating protein MSTP9 (Brain rescue factor 1) (BRF-1) (Hepatocyte growth factor-like protein homolog)
Gene Name	MST1P9 D1F15S1A MSTP9
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q2TV78
Alternative Names	
Background	

