



MFS2B rabbit pAb

Cat No.:ES19848

For research use only

Overview

Product Name	MFS2B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MFS2B AA range: 145-195
Specificity	This antibody detects endogenous levels of human MFS2B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Major facilitator superfamily domain-containing protein 2B
Gene Name	MFSD2B
Cellular localization	Cell membrane ; Multi-pass membrane protein . Localizes to the cell membrane and intracellular membranes. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	388931
Human Swiss-Prot Number	A6NFX1
Alternative Names	
Background	

