



YPEL3 rabbit pAb

Cat No.:ES19854

For research use only

Overview

Product Name	YPEL3 rabbit pAb
Host species	Rabbit
Applications	IHC
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	IHC-p 1: 50-200
Immunogen	Synthesized peptide derived from human YPEL3. AA range:100-119
Specificity	This antibody detects endogenous levels of human YPEL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Protein yippee-like 3
Gene Name	YPEL3 FKSG5
Cellular localization	Nucleus, nucleolus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	83719
Human Swiss-Prot Number	P61236
Alternative Names	
Background	

